




products for micro surgery super micro-surgery



Disposable Micro Vascular Clamps

KEISEI VASCULAR CLAMP												
ordering code number	S11Y	A11Y	S11B	A11B	S22Y	A22Y	S22B	A22B	S33Y	A33Y	S33B	A33B
color Yellow / Blue	Y	Y	B	B	Y	Y	B	B	Y	Y	B	B
single / double	s	d	s	d	s	d	s	d	s	d	s	d
Venous / Artery	V		V/A		V		V/A		V		V/A	
clamp pressure (gf)	8 - 15		13.5 - 20		15 - 25		20 - 23.5		30 - 50		45 - 60	
applicable vessel diameter (mm)	0.4 - 1.0				0.6 - 1.5				1.0 - 2.25			
	A	12.6mm			15.3mm			17.3mm				
	B	4.5mm			5mm			5mm				
	C	2.6mm			3mm			3.5mm				
	D	16mm			20mm			21mm				

Disposable Background material

25mm x 50mm
12pcs/ pack



CB1



CB2



CB3



CB4

Disposable Ruler

MI 1847-D 25pcs/ pack



KEISEI MEDICAL

Indian Distributors:

Genetic

Genetic Instruments Pvt. Ltd.,

Plot No.2, Sree Krishna Nagar 3rd Street, Kolathur, Chennai - 600 099.

Tel: +91-44-2651 1506 / 2651 1706 / 2651 1806 E-mail: genetic@geneticinstruments.com

Microsurgery Instruments



① **TF-51C 11cm**
Micro Forceps


② **MF-05 11cm**
Micro Forceps

③ **TF-5 11cm**
Micro Forceps

④ **VF-5 11cm**
Vascular Dilator

⑤ **AMF-1 12cm**
Micro Adson Forceps With hook





Sterile Packing

PRODUCT

Mini - Red	Maxi - Red
Mini - Blue	Maxi - Blue
Mini - Yellow	Maxi Yellow
Mini - White	Maxi - White

Colour Code

Red - Arteries	Blue - Veins
Yellow - Ureter	White - Nerves and Tendons

Identi-Loops are medical grade silicone rubber ribbons useful in surgical situations requiring retraction, identification or occlusion of arteries, veins, nerves, ureters or tendons.

Design

Identi-Loops are oval in cross-section allowing for even pressure distribution, thus minimizing damage to tissue.

Maximum Safety

Identi-Loops are made of Silclear Silicone rubber, a smooth, soft, bio-compatible material. The Identi-Loops are non-thrombogenic, non-leeachable, non-toxic and resist oxidation. They do not support bacterial growth. Identi-Loops are radiopaque and lot traceable to ensure patient safety. Non-absorbent and non-adherent to tissue and to themselves, they are pliable and will not dry to a hard edge, thus enabling quick and safe release.

Identi-Loops are manufactured in a clean room environment, meeting the US Pharmacopoeia XX Class VI, Good Manufacturing Practice, PDA and other international standards.

Packaging

Identi-Loops are packaged sterile, Box of 10 Pouches, (2 Loops per Pouch)

Genetic

GENETIC INSTRUMENTS PVT LTD

Plot No.2, Sri Krishna Nagar, 3rd Street, Kolathur, Chennai - 600 099.

Ph: +91-44-26511506, 26511706, 26511806

e-mail: genetic@geneticinstruments.com